



I M P L A N T S

I N T E R N A T I O N A L

B I O L O G I C S

## CRANIO-FIX™

### A MODERN ACRYLIC BONE CEMENT

#### PROPERTIES:

- Standard Viscosity
- Radio Opacity
- Lower proportions of monomer released in situ
- Improved tactile feel
- Reduced curing time
- Rapid final hardening phase, reducing the risk of prosthetic micro-movement
- Lower polymerisation temperature
- Mechanical properties far exceed the ISO 5833 standard

#### INDICATIONS:

- Our Acrylic Bone Cement preparation can be used to fix a metallic implant to living bone and or to repair defects in the cranium caused by either trauma or other pathologies

#### PACKAGING:


- Liquid contained within glass ampoule - 14.4g.
- Powder contained within double Tyvec® pouch - 40g.

#### STORAGE:

- This product is delivered sterile, ready for use in theatre and must be stored in a clean, dark and dry location at an ambient temperature of between 5°C and 25°C

#### ACCESSORIES:

- Mixing bowl and spatula available on demand



#### CATALOGUE #:

- Cranio-Fix™ Cement BC0003
- Mixing bowl BC0012
- Spatula BC0013