



# IMPLANTS

INTERNATIONAL

ARTHROPLASTY



**Wagner Moduloc™**  
Fixation Bolt

**Wagner Moduloc™**  
Prosthesis Body

**Wagner Moduloc™**  
Stem



**Wagner Moduloc™**  
Modular Revision  
Hip Stem

**WAGNER MODULOC™**  
MODULAR REVISION HIP STEM



# Wagner Moduloc™ Revision Hip Stem

## Wagner Moduloc™ Revision Hip Prosthesis Fixation Bolt

Size Length(L1 mm) Cat#

7.00	29	HS1415
8.00	35	HS1416
9.00	45	HS1417
10.00	55	HS1418

## Wagner Moduloc™ Revision Hip Prosthesis Body with Tendon Holes

Size Length (L2mm) Cat#

X-Short	42	HS1419
Short	48	HS1420
Medium	58	HS1421
Long	68	HS1422

## Wagner Moduloc™ Revision Hip Prosthesis Body Standard

Size Length (L3mm) Cat#


X-Short	42	HS1423
Short	48	HS1424
Medium	58	HS1425
Long	68	HS1426

## Wagner Moduloc™ Revision Hip Prosthesis Stem

Stem Size (Dmm) Stem Length (L4mm) Taper Cat#

13	142	12/14	HS1427
13	172	12/14	HS1428
13	212	12/14	HS1429
14	142	12/14	HS1430
14	172	12/14	HS1431
14	212	12/14	HS1432
15	142	12/14	HS1433
15	172	12/14	HS1434
15	212	12/14	HS1435
16	142	12/14	HS1436
16	172	12/14	HS1437
16	212	12/14	HS1438
17	142	12/14	HS1439
17	172	12/14	HS1440
17	212	12/14	HS1441
18	142	12/14	HS1442
18	172	12/14	HS1443
18	212	12/14	HS1444
19	142	12/14	HS1445
19	172	12/14	HS1446
19	212	12/14	HS1447
20	142	12/14	HS1448
20	172	12/14	HS1449
20	212	12/14	HS1450

Material: Titanium Alloy to ISO5832-3



September 2013

Manufactured in the UK

DISTRIBUTED BY

Implants International Ltd  
 71 Jay Avenue, Teesside Industrial Estate  
 Thornaby-on-Tees, TS17 9LZ, United Kingdom  
 tel: +44(0)1642 769080 fax: +44(0)1642 765848  
 e-mail: enquiries@implantsinternational.com  
 web: www.implantsinternational.com