



Low Profile Anterior Cervical Plate

Level 2 - 38mm




1.8 / 1.9mm Thick


Plates Sizes Available

- Level 1 - 6 Sizes from 20mm - 30mm in 2mm increments
- Level 2 - 6 Sizes from 38mm - 48mm in 2mm increments
- Level 3 - 5 Sizes from 54mm - 70mm in 4mm increments
- Level 4 - 5 Sizes from 84mm - 100mm in 2mm increments

Level 2 - 38mm



R400mm Lordosis



R40mm Axial



**Implants
International**

71 Jay Avenue Thornaby-on-Tees
North Riding of Yorkshire
TS17 9LZ

United Kingdom

tel: 01642 769080

email: enquiries@implantsinternational.com

web: www.implantsinternational.com