


Modular Foot Implant Systems



Contents

Sterile product and the dedicated instrument sets are supplied in their own cases to improve CSSD turnaround time and reduce the cost of sterilizing sets not used in surgery

Each Module comprises a Sterile Store and a Dedicated Instrument Set

- Module 1** TARGET™ Type 1 & 2 Cannulated Compression Screws & Harpoon™
- S³ Screws
 - Percutaneous Fracture / Fusion Harpoon™
 - TARGET™ Screws - Optional Sterile Disposable Surgery Pack

- Module 2** Combined Serrated Staples & Locking Compression Plates
- MEGA™ Compression Locking Plates
 - MEGA™ Optional Sterile Disposable Surgery Pack
 - Compression Staple
 - Akin Staple

- Module 3** Anatomic MTP Fusion Plates: Non-locking
- D-Tach™ Snap-off Screws
 - D-Tach™ Weil Snap-off Screws
 - Standard Plates – Non-locking
 - Revision Plates – Non-locking
 - MTP Fusion – Optional Sterile Disposable Surgery Pack

- Module 4** Locking Plates and Screws: Fore, Mid & Hind Foot indications
- Small Extremity Plates and Screw Systems

- Module 5** HV Xpanding Spring-Clip™: IM implant for Hallux Valgus Deformity Correction

- Module 6** XpandComp™
- Hind Foot (HF) Fusion Device
 - Tibial Talar Calcaneal (TTC) Fusion Rod

- Module 7** Minimally Invasive System
- Burrs and Cutters
 - Instruments

- Module 8** Wire Compression & Distraction Clamp / Pad Distractor Clamp

Appendices


Osteotomy / Scarf Instrument Set Typical Module 1



Dedicated Outer Storage & Carry Case




System Instrument Tray




TARGET™ Drill and K-Wire Guide on inside lid

TARGET™
Cannulated Compression Screw Sterile Stores




Dedicated Outer Storage and Carry Case




Compartmentalized Inner Tray

TARGET™ S³ Compression Screw Osteotomy Set



Instrument Set # IN9854 DI1001 & DI1002

TARGET™ S³ Osteotomy Set



Cannulated self-tapping and self-drilling Compression Screws
with anti-backout feature and opposing buttress threads


TARGET™ 1 S³ Ø3mm x 10mm to 40mm in 2mm increments plus an extra-long length of 46mm.
Used mainly for fore-foot indications.

TARGET™ 2 S³ Ø4mm x 12mm to 34mm long in 2mm increments.
Used mainly in mid-foot indications but can be used for forefoot in certain circumstances
when greater strength and stability are required.

Both screws have Ø1.1mm cannulations to work with 1mm x 80mm long K-Wires.

- We offer 3 lengths of 2 1 Cannulated Drills and Counter-sinks for Type 1 Screws marked A, B, and C to suit various screw lengths. For user details, please refer to the attached Drill Chart: [MD163 Rev.0]
- We offer 3 lengths of 2 1 Cannulated Drills and Counter-sinks for Type 2 screws marked D, E, and F to suit various screw lengths. For user details, please refer to the attached Drill Chart: [MD163 Rev.0]
- We also offer two sets of sterile end-cutting Cannulated Drills and separate Cannulated Chamfer Cutters to suit TARGET™ 1 & TARGET™ 2

Our set includes a Chevron Jig, designed to offer accurate 90° (<14° deformity corrections) and 60° (>14°) osteotomies, supplied with a dedicated saw blade and Ø1.5mm K-Wires.




Included in the set is an Osteotomy Bone Clamp to help stabilize the osteotomy whilst the site is drilled and stabilized with screws and or staples.

Harpoon™ Percutaneous Fixation and Fusion Device

Suitable for a miscellany of indications including, but not limited to:

- Hammer Toes
- Claw Toes
- Mallet Toes
- Metatarsal Fractures
- ✓ Lisfranc fractures
- ✓ Midfoot Arthrodesis
- ✓ Forefoot Arthrodesis



Ø2.5mm x 150mm long with 5mm, 10mm, 20mm offsets and dual trocar tips

Trocar tipped



Standard Compression Nut
TS8125



Flanged Compression Nut
For use with cerclage wire
TS8124



- No pre-drilling required
- Dual trocar ended: suitable for retrograde entry, followed by antegrade advance, placement, and compression
- Optimal fracture site compression with a choice of compression fittings
- Non-theatre retrieval of device

Locking Compression Plates

□ MEGA™ Compression Staple Locking Plates: DUO & QUAD

These are Low-Profile (1.5mm thick) Compression Plates designed for use with our Ø3mm self-drilling and self-tapping taper lock screws. MEGA staples may be combined together to span greater lengths when required.

The screw hole centres on these staples are offset to minimize stress shielding of screws deployed. Screw heads are designed to taper-lock into each screw hole. These devices are suitable for midfoot or hindfoot indications. For further details, please refer to the attached user chart [MD162 Rev.0].

A drill jig is provided in the set to facilitate a 1.9mm pilot hole to receive the Ø3mm screws. Use our compression staple adjustable drill jig as a size gauge.

Duo



15mm, 20mm & 25mm
Hole – IA centres
Ref. MD162 Rev.0

Quad



□ MEGA™ Optional Sterile Disposable Surgery Pack




Sterile-wrapped and housed in protective blister packs

Compression Serrated Staples


Serrated rectangular-sectioned Compression Staple

Ceramic-bead peened to stress-relieve and prevent recoil, post distraction. These can be used in mid-foot indications. For details of IA and leg lengths, please refer to MD161 Rev. 0. Pre-drill with the 2.2mm pilot drill and adjustable drill jig provided in the instrument set.



Serrated square-sectioned Akin Staple

Sized at 1 x 1 x 9mm, used mainly for forefoot indications where you predrill with the Ø1.5mm K-wire and dedicated drill jig provided in the instrument set.




MTP Fusion non-locking set

We offer 3 lengths each of Ultra-Low Profile Anatomic Plates designed for dorsal placement. These are non-locking plates designed for use with our D-Tach™ screws.

We also offer 4 matching diameters (Ø14mm, Ø16mm, Ø18mm, and Ø20mm) of Cannulated Hemispherical Cutters designed to create close-fitting hemispherical domes for the metatarsal head and base of Phalanx.

Fixation is via the use of Ø2.7mm, Ø3mm, or Ø2mm Weil D-Tach™ snap-off screws. The Ø3mm screws are designed for rescue or more demanding situations. These screws require pre-drilling in hard bone with Ø1.5mm or Ø1.9mm dedicated drills. When a pilot hole is deemed necessary, our drill jigs must be used. These are designed to locate within the screw holes in the plate. We strongly recommend that our drill jigs be used to pre-drill the bone, as free-hand drilling with Ø1.9 mm drills will over-drill holes.



If the drill jig is not used, then free-hand drill with the Ø1.9mm drill provided. It is also possible to gently power-drive our D-Tach™ screws without pre-drilling the bone. However, if this method is preferred, it is prudent to leave the final torqueing by hand after the detachable drives are gently detached, whilst holding the plate to prevent de-stabilizing the construct.

We also offer a graduated plate bender for use with our Standard anatomic plates.

Standard Plates may be formed to suit the patient's metatarsal. However, these Titanium plates MUST NOT be overbent and corrected, as that will weaken these plates.



Max Plates are supplied **pre-formed**. These plates must not be manipulated or re-formed.



D-Tach™ screws may be advanced under power 95% of the way in, but final tightening MUST be undertaken using our dedicated drivers.

The Plate needs to be held securely in situ before detaching the snap-off drive shafts of our D-Tach™ screws, once all of them have been tightened. Holding the plate will minimize the risk of destabilizing the screw during the detaching of the drive tab.

HV (Hallux Valgus) Xpanding Spring-Clip™

Intramedullary alignment and fixation device for treating Hallux Valgus Deformities

Pre-formed high tensile stainless-steel expanding spring clips offered with a choice of 2mm, 4mm and 6mm offsets. Using a standard medial approach and the dedicated anatomic resection guide to ensure optimum anatomic and linear correction of the metatarsal head.


HV Xpanding Spring-Clip™

2mm Offset




TS8375

4mm Offset




TS8376

6mm Offset




TS8377

Locking Screw




CAT #	Size	CAT #	Size
TS8378	Ø4 x 17mm	TS8381	Ø4 x 20mm
TS8379	Ø4 x 18mm	TS8382	Ø4 x 21mm
TS8380	Ø4 x 19mm	TS8383	Ø4 x 22mm

Left Resection



IN3025


Right Resection Jig



IN3025

XPANDCOMP™ Systems

XPANDCOMP™ Hindfoot Fusion Device




“Ø” Refers to pre-expansion



Available lengths: 60mm, 70mm, 80mm, 90mm, 100mm, 110mm

XPANDCOMP™ Tibial Talar Calcaneal Fusion Device



“Ø” Refers to pre-expansion



Available lengths: 120mm, 140mm, 160mm, 180mm, 200mm, 220mm

Cold Formed, Low-Profile Titanium Non-Rigid Locking Plates Fore, Mid & Hind Foot Indications

Customized Sets to Suit Usage at each Hospital

Plates with Colour Coded Screw Holes to Facilitate Screw Ø

Selection. Refer to Pages 14 to 17 for Catalogue Numbers

**MTP1 Fusion 8°
Dorsal Approach**




LH / RH
S, M, L

**Metatarsal Basal Osteotomy
Chevron Type I (Small bones)**



LH / RH
S, M, L

**Metatarsal Basal Osteotomy
Chevron Type II (Large bones)**




LH / RH
S, M, L

**Metatarsal Osteotomy
Wedge for Chevron Type I & II**



4, 5, 6,
7mm

**Metatarsal Osteotomy
Ludloff**




S, M, L

**Metatarsal MIS
Hallux Valgus Correction**




LH / RH
S, L

**TMT1 Fusion
Medial Approach**




LH / RH
S, L

**TMT1 Fusion
Plantar Approach**



LH / RH

**L Plate
(8H Shown)**



LH / RH
4H, 6H, 8H

**T Plate
(6H Shown)**



4H, 6H

X Plate




S, M, L

**Straight Plate
(6H Shown)**




4H, 6H, 8H

Cuboid Plate




LH / RH

Grid Type Plate




S, M, L

Navicular Plate




LH / RH

Talus Plate




LH / RH

Calcaneal
Offset




3, 5, 7,
9mm

Calcaneal
Type I
(LH, S Shown)



LH / RH
S, L


Calcaneal
Type II
(RH, S Shown)



LH / RH
S, L


Screws

Locking Screw
Ø2.4 mm, T10




10 – 30mm (2mm increments)

Locking Screw
Ø3.2 mm, T15




10 – 40mm (2mm increments)

Cortical Screw
Ø2.4 mm, T10



10 – 30mm (2mm increments)

Cortical Screw
Ø3.2 mm, T15



10 – 40mm (2mm increments)

CAT #	Item	Type	Size	Material
FL2001	MTP1 FUSION 8° PLATE	DORSAL APPROACH	LEFT SHORT	Ti. 6/4 / ISO:5832-3
FL2002	MTP1 FUSION 8° PLATE	DORSAL APPROACH	LEFT MEDIUM	Ti. 6/4 / ISO:5832-3
FL2003	MTP1 FUSION 8° PLATE	DORSAL APPROACH	LEFT LONG	Ti. 6/4 / ISO:5832-3
FL2004	MTP1 FUSION 8° PLATE	DORSAL APPROACH	RIGHT SHORT	Ti. 6/4 / ISO:5832-3
FL2005	MTP1 FUSION 8° PLATE	DORSAL APPROACH	RIGHT MEDIUM	Ti. 6/4 / ISO:5832-3
FL2006	MTP1 FUSION 8° PLATE	DORSAL APPROACH	RIGHT LONG	Ti. 6/4 / ISO:5832-3
FL2007	CHEVRON TYPE I PLATE	METATARSAL BASAL OSTEOTOMY	LEFT SHORT	Ti. 6/4 / ISO:5832-3
FL2008	CHEVRON TYPE I PLATE	METATARSAL BASAL OSTEOTOMY	LEFT MEDIUM	Ti. 6/4 / ISO:5832-3
FL2009	CHEVRON TYPE I PLATE	METATARSAL BASAL OSTEOTOMY	LEFT LONG	Ti. 6/4 / ISO:5832-3
FL2010	CHEVRON TYPE I PLATE	METATARSAL BASAL OSTEOTOMY	RIGHT SHORT	Ti. 6/4 / ISO:5832-3
FL2011	CHEVRON TYPE I PLATE	METATARSAL BASAL OSTEOTOMY	RIGHT MEDIUM	Ti. 6/4 / ISO:5832-3
FL2012	CHEVRON TYPE I PLATE	METATARSAL BASAL OSTEOTOMY	RIGHT LONG	Ti. 6/4 / ISO:5832-3
FL2013	CHEVRON TYPE II PLATE	METATARSAL BASAL OSTEOTOMY	LEFT SHORT	Ti. 6/4 / ISO:5832-3
FL2014	CHEVRON TYPE II PLATE	METATARSAL BASAL OSTEOTOMY	LEFT MEDIUM	Ti. 6/4 / ISO:5832-3
FL2015	CHEVRON TYPE II PLATE	METATARSAL BASAL OSTEOTOMY	LEFT LONG	Ti. 6/4 / ISO:5832-3
FL2016	CHEVRON TYPE II PLATE	METATARSAL BASAL OSTEOTOMY	RIGHT SHORT	Ti. 6/4 / ISO:5832-3
FL2017	CHEVRON TYPE II PLATE	METATARSAL BASAL OSTEOTOMY	RIGHT MEDIUM	Ti. 6/4 / ISO:5832-3
FL2018	CHEVRON TYPE II PLATE	METATARSAL BASAL OSTEOTOMY	RIGHT LONG	Ti. 6/4 / ISO:5832-3
FL2019	WEDGE FOR CHEVRON TYPE I & II	METATARSAL BASAL OSTEOTOMY	4MM	Ti. 6/4 / ISO:5832-3
FL2020	WEDGE FOR CHEVRON TYPE I & II	METATARSAL BASAL OSTEOTOMY	5MM	Ti. 6/4 / ISO:5832-3
FL2021	WEDGE FOR CHEVRON TYPE I & II	METATARSAL BASAL OSTEOTOMY	6MM	Ti. 6/4 / ISO:5832-3
FL2022	WEDGE FOR CHEVRON TYPE I & II	METATARSAL BASAL OSTEOTOMY	7MM	Ti. 6/4 / ISO:5832-3
FL2023	LUDLOFF PLATE	METATARSAL OSTEOTOMY - MEDIAL	SHORT	Ti. 6/4 / ISO:5832-3
FL2024	LUDLOFF PLATE	METATARSAL OSTEOTOMY - MEDIAL	MEDIUM	Ti. 6/4 / ISO:5832-3
FL2025	LUDLOFF PLATE	METATARSAL OSTEOTOMY - MEDIAL	LONG	Ti. 6/4 / ISO:5832-3
FL2026	HALLUX VALGUS CORRECTION PLATE	METATARSAL MIS MEDIAL APPROACH	LEFT SHORT	Ti. 6/4 / ISO:5832-3
FL2027	HALLUX VALGUS CORRECTION PLATE	METATARSAL MIS MEDIAL APPROACH	LEFT LONG	Ti. 6/4 / ISO:5832-3
FL2028	HALLUX VALGUS CORRECTION PLATE	METATARSAL MIS MEDIAL APPROACH	RIGHT SHORT	Ti. 6/4 / ISO:5832-3
FL2029	HALLUX VALGUS CORRECTION PLATE	METATARSAL MIS MEDIAL APPROACH	RIGHT LONG	Ti. 6/4 / ISO:5832-3
FL2030	TMT1 FUSION PLATE	MEDIAL APPROACH	LEFT SHORT	Ti. 6/4 / ISO:5832-3
FL2031	TMT1 FUSION PLATE	MEDIAL APPROACH	LEFT LONG	Ti. 6/4 / ISO:5832-3
FL2032	TMT1 FUSION PLATE	MEDIAL APPROACH	RIGHT SHORT	Ti. 6/4 / ISO:5832-3
FL2033	TMT1 FUSION PLATE	MEDIAL APPROACH	RIGHT LONG	Ti. 6/4 / ISO:5832-3

FL2034	TMT1 FUSION PLATE	PLANTAR APPROACH	LEFT	Ti. 6/4 / ISO:5832-3
FL2035	TMT1 FUSION PLATE	PLANTAR APPROACH	RIGHT	Ti. 6/4 / ISO:5832-3
FL2036	L PLATE	DORSAL APPROACH	LEFT 4H	Ti. 6/4 / ISO:5832-3
FL2037	L PLATE	DORSAL APPROACH	LEFT 6H	Ti. 6/4 / ISO:5832-3
FL2038	L PLATE	DORSAL APPROACH	LEFT 8H	Ti. 6/4 / ISO:5832-3
FL2039	L PLATE	DORSAL APPROACH	RIGHT 4H	Ti. 6/4 / ISO:5832-3
FL2040	L PLATE	DORSAL APPROACH	RIGHT 6H	Ti. 6/4 / ISO:5832-3
FL2041	L PLATE	DORSAL APPROACH	RIGHT 8H	Ti. 6/4 / ISO:5832-3
FL2042	T PLATE	DORSAL APPROACH	4H	Ti. 6/4 / ISO:5832-3
FL2043	T PLATE	DORSAL APPROACH	6H	Ti. 6/4 / ISO:5832-3
FL2044	X PLATE	DORSAL APPROACH	SHORT	Ti. 6/4 / ISO:5832-3
FL2045	X PLATE	DORSAL APPROACH	MEDIUM	Ti. 6/4 / ISO:5832-3
FL2046	X PLATE	DORSAL APPROACH	LONG	Ti. 6/4 / ISO:5832-3
FL2047	STRAIGHT PLATE	DORSAL APPROACH	4H	Ti. 6/4 / ISO:5832-3
FL2048	STRAIGHT PLATE	DORSAL APPROACH	6H	Ti. 6/4 / ISO:5832-3
FL2049	STRAIGHT PLATE	DORSAL APPROACH	8H	Ti. 6/4 / ISO:5832-3
FL2050	CUBOID PLATE	DORSAL APPROACH	LEFT	Ti. 6/4 / ISO:5832-3
FL2051	CUBOID PLATE	DORSAL APPROACH	RIGHT	Ti. 6/4 / ISO:5832-3
FL2052	GRID TYPE PLATE	DORSAL APPROACH	SHORT	Ti. 6/4 / ISO:5832-3
FL2053	GRID TYPE PLATE	DORSAL APPROACH	MEDIUM	Ti. 6/4 / ISO:5832-3
FL2054	GRID TYPE PLATE	DORSAL APPROACH	LONG	Ti. 6/4 / ISO:5832-3
FL2055	NAVICULAR PLATE	DORSAL APPROACH	LEFT	Ti. 6/4 / ISO:5832-3
FL2056	NAVICULAR PLATE	DORSAL APPROACH	RIGHT	Ti. 6/4 / ISO:5832-3
FL2057	TALUS PLATE	DORSAL APPROACH	LEFT	Ti. 6/4 / ISO:5832-3
FL2058	TALUS PLATE	DORSAL APPROACH	RIGHT	Ti. 6/4 / ISO:5832-3
FL2059	CALCANEAL OFFSET	LATERAL APPROACH	3MM	Ti. 6/4 / ISO:5832-3
FL2060	CALCANEAL OFFSET	LATERAL APPROACH	5MM	Ti. 6/4 / ISO:5832-3
FL2061	CALCANEAL OFFSET	LATERAL APPROACH	7MM	Ti. 6/4 / ISO:5832-3
FL2062	CALCANEAL OFFSET	LATERAL APPROACH	9MM	Ti. 6/4 / ISO:5832-3
FL2063	CALCANEAL TYPE I PLATE	LATERAL APPROACH	LEFT SHORT	Ti. 6/4 / ISO:5832-3
FL2064	CALCANEAL TYPE I PLATE	LATERAL APPROACH	LEFT LONG	Ti. 6/4 / ISO:5832-3
FL2065	CALCANEAL TYPE I PLATE	LATERAL APPROACH	RIGHT SHORT	Ti. 6/4 / ISO:5832-3
FL2066	CALCANEAL TYPE I PLATE	LATERAL APPROACH	RIGHT LONG	Ti. 6/4 / ISO:5832-3
FL2067	CALCANEAL TYPE II PLATE	LATERAL APPROACH	LEFT SHORT	Ti. 6/4 / ISO:5832-3
FL2068	CALCANEAL TYPE II PLATE	LATERAL APPROACH	LEFT LONG	Ti. 6/4 / ISO:5832-3
FL2069	CALCANEAL TYPE II PLATE	LATERAL APPROACH	RIGHT SHORT	Ti. 6/4 / ISO:5832-3
FL2070	CALCANEAL TYPE II PLATE	LATERAL APPROACH	RIGHT LONG	Ti. 6/4 / ISO:5832-3

CAT #	Item	Size
FL2101	CORTICAL SCREW	Ø2.4 x 10MM
FL2102	CORTICAL SCREW	Ø2.4 x 12MM
FL2103	CORTICAL SCREW	Ø2.4 x 14MM
FL2104	CORTICAL SCREW	Ø2.4 x 16MM
FL2105	CORTICAL SCREW	Ø2.4 x 18MM
FL2106	CORTICAL SCREW	Ø2.4 x 20MM
FL2107	CORTICAL SCREW	Ø2.4 x 22MM
FL2108	CORTICAL SCREW	Ø2.4 x 24MM
FL2109	CORTICAL SCREW	Ø2.4 x 26MM
FL2110	CORTICAL SCREW	Ø2.4 x 28MM
FL2111	CORTICAL SCREW	Ø2.4 x 30MM

FL2201	CORTICAL SCREW	Ø3.2 x 10MM
FL2202	CORTICAL SCREW	Ø3.2 x 12MM
FL2203	CORTICAL SCREW	Ø3.2 x 14MM
FL2204	CORTICAL SCREW	Ø3.2 x 16MM
FL2205	CORTICAL SCREW	Ø3.2 x 18MM
FL2206	CORTICAL SCREW	Ø3.2 x 20MM
FL2207	CORTICAL SCREW	Ø3.2 x 22MM
FL2208	CORTICAL SCREW	Ø3.2 x 24MM
FL2209	CORTICAL SCREW	Ø3.2 x 26MM
FL2210	CORTICAL SCREW	Ø3.2 x 28MM
FL2211	CORTICAL SCREW	Ø3.2 x 30MM
FL2212	CORTICAL SCREW	Ø3.2 x 32MM
FL2213	CORTICAL SCREW	Ø3.2 x 34MM
FL2214	CORTICAL SCREW	Ø3.2 x 36MM
FL2215	CORTICAL SCREW	Ø3.2 x 38MM
FL2216	CORTICAL SCREW	Ø3.2 x 40MM

CAT #	Item	Size
FL2151	LOCKING CORTICAL SCREW	Ø2.4 x 10MM
FL2152	LOCKING CORTICAL SCREW	Ø2.4 x 12MM
FL2153	LOCKING CORTICAL SCREW	Ø2.4 x 14MM
FL2154	LOCKING CORTICAL SCREW	Ø2.4 x 16MM
FL2155	LOCKING CORTICAL SCREW	Ø2.4 x 18MM
FL2156	LOCKING CORTICAL SCREW	Ø2.4 x 20MM
FL2157	LOCKING CORTICAL SCREW	Ø2.4 x 22MM
FL2158	LOCKING CORTICAL SCREW	Ø2.4 x 24MM
FL2159	LOCKING CORTICAL SCREW	Ø2.4 x 26MM
FL2160	LOCKING CORTICAL SCREW	Ø2.4 x 28MM
FL2161	LOCKING CORTICAL SCREW	Ø2.4 x 30MM

FL2251	LOCKING CORTICAL SCREW	Ø3.2 x 10MM
FL2252	LOCKING CORTICAL SCREW	Ø3.2 x 12MM
FL2253	LOCKING CORTICAL SCREW	Ø3.2 x 14MM
FL2254	LOCKING CORTICAL SCREW	Ø3.2 x 16MM
FL2255	LOCKING CORTICAL SCREW	Ø3.2 x 18MM
FL2256	LOCKING CORTICAL SCREW	Ø3.2 x 20MM
FL2257	LOCKING CORTICAL SCREW	Ø3.2 x 22MM
FL2258	LOCKING CORTICAL SCREW	Ø3.2 x 24MM
FL2259	LOCKING CORTICAL SCREW	Ø3.2 x 26MM
FL2260	LOCKING CORTICAL SCREW	Ø3.2 x 28MM
FL2261	LOCKING CORTICAL SCREW	Ø3.2 x 30MM
FL2262	LOCKING CORTICAL SCREW	Ø3.2 x 32MM
FL2263	LOCKING CORTICAL SCREW	Ø3.2 x 34MM
FL2264	LOCKING CORTICAL SCREW	Ø3.2 x 36MM
FL2265	LOCKING CORTICAL SCREW	Ø3.2 x 38MM
FL2266	LOCKING CORTICAL SCREW	Ø3.2 x 40MM

R. Small Extremity Plates and Screw Systems

Ref No	EAN13 Barcode	Description
Small Extremity Plates – Forearm (Curved, Titanium)		
152.0021.1006	8682241800971	Plate, Small Forearm, Curved, Ti. 6 Holes
152.0021.1007	8682241800988	Plate, Small Forearm, Curved, Ti. 7 Holes
152.0021.1008	8682241800995	Plate, Small Forearm, Curved, Ti. 8 Holes
152.0021.1009	8682241801008	Plate, Small Forearm, Curved, Ti. 9 Holes
152.0021.1010	8682241801015	Plate, Small Forearm, Curved, Ti. 10 Holes
152.0021.1011	8682241801022	Plate, Small Forearm, Curved, Ti. 11 Holes
152.0021.1012	8682241801039	Plate, Small Forearm, Curved, Ti. 12 Holes
Small Extremity Plates – Forearm (Straight, Titanium)		
152.0022.1006	8682241801046	Plate, Small Forearm, Straight, Ti. 6 Holes
152.0022.1007	8682241801053	Plate, Small Forearm, Straight, Ti. 7 Holes
152.0022.1008	8682241801060	Plate, Small Forearm, Straight, Ti. 8 Holes
152.0022.1009	8682241801077	Plate, Small Forearm, Straight, Ti. 9 Holes
152.0022.1010	8682241801084	Plate, Small Forearm, Straight, Ti. 10 Holes
152.0022.1011	8682241801091	Plate, Small Forearm, Straight, Ti. 11 Holes
152.0022.1012	8682241801107	Plate, Small Forearm, Straight, Ti. 12 Holes
Distal Radius Volar Plates – Small (Left/Right)		
152.0025.1103	8682241801114	Plate, Distal Radius Volar, Ti. 3 Holes, Small, Left
152.0025.1104	8682241801121	Plate, Distal Radius Volar, Ti. 4 Holes, Small, Left
152.0025.1105	8682241801138	Plate, Distal Radius Volar, Ti. 5 Holes, Small, Left
152.0025.1203	8682241801145	Plate, Distal Radius Volar, Ti. 3 Holes, Small, Right
152.0025.1204	8682241801152	Plate, Distal Radius Volar, Ti. 4 Holes, Small, Right
152.0025.1205	8682241801169	Plate, Distal Radius Volar, Ti. 5 Holes, Small, Right
152.0026.1103	8682241801176	Plate, Distal Radius Volar, Ti. 3 Holes, Left
152.0026.1104	8682241801183	Plate, Distal Radius Volar, Ti. 4 Holes, Left
152.0026.1105	8682241801190	Plate, Distal Radius Volar, Ti. 5 Holes, Left
152.0026.1203	8682241801206	Plate, Distal Radius Volar, Ti. 3 Holes, Right
152.0026.1204	8682241801213	Plate, Distal Radius Volar, Ti. 4 Holes, Right
152.0026.1205	8682241801220	Plate, Distal Radius Volar, Ti. 5 Holes, Right
Distal Radius Dorsal Plates (Left & Right)		
152.0027.1103	8682241801237	Plate, Distal Radius Dorsal, Left, Ti., 3 D
152.0027.1104	8682241801244	Plate, Distal Radius Dorsal, Left, Ti., 4 D
152.0027.1105	8682241801251	Plate, Distal Radius Dorsal, Left, Ti., 5 D
152.0027.1203	8682241801268	Plate, Distal Radius Dorsal, Right, Ti., 3 D
152.0027.1204	8682241801275	Plate, Distal Radius Dorsal, Right, Ti., 4 D
152.0027.1205	8682241801282	Plate, Distal Radius Dorsal, Right, Ti., 5 D
Distal Hook Plates		
152.0028.1003	8682241801299	Plate, Distal Type Hook, Ti. 3 Holes
152.0028.1004	8682241801305	Plate, Distal Type Hook, Ti. 4 Holes
152.0028.1005	8682241801312	Plate, Distal Type Hook, Ti. 5 Holes
Ulnar Tip Plates (Left & Right)		
152.0029.1103	8682241801329	Plate, Ulnar Tip, Left, Ti., 3 D
152.0029.1104	8682241801336	Plate, Ulnar Tip, Left, Ti., 4 D

152.0029.1105	8682241801343	Plate, Ulnar Tip, Left, Ti., 5 D
152.0029.1203	8682241801350	Plate, Ulnar Tip, Right, Ti., 3 D
152.0029.1204	8682241801367	Plate, Ulnar Tip, Right, Ti., 4 D
152.0029.1205	8682241801374	Plate, Ulnar Tip, Right, Ti., 5 D

Hand Locking Plates

152.0031.1528	8682241840045	Locking Plate, Hand, straight, 10 holes
152.0031.1529	8682241840052	Locking Plate, Hand, straight compression, 4 holes
152.0031.1530	8682241840069	Locking Plate, Hand, straight compression, 6 holes
152.0031.1531	8682241840076	Locking Plate, Hand, T type, 8+3 holes
152.0031.1532	8682241840083	Locking Plate, Hand, T type, 8+4 holes
152.0031.1533	8682241840090	Locking Plate, Hand, Y type, 8+3 holes
152.0031.1534	8682241840106	Locking Plate, Hand, L type, 4+2 holes, Left angled
152.0031.1535	8682241840113	Locking Plate, Hand, L type, 8+2 holes, Left angled
152.0031.1536	8682241840120	Locking Plate, Hand, L type, 4+2 holes, Right angled
152.0031.1537	8682241840137	Locking Plate, Hand, L type, 8+2 holes, Right angled
152.0031.1538	8682241840144	Locking Plate, Hand, Condylar Y type, 6+2 holes
152.0031.1542	8682241840151	Locking Plate, Hand, Scaphoid, 6 holes
152.0031.1543	8682241840168	Locking Plate, Hand, Scaphoid, 6 holes
152.0031.1544	8682241840175	Locking Plate, Hand, Rotation correction, 4+2 holes
152.0031.1545	8682241840182	Locking Plate, Hand, Rotation correction, 4+4 holes
152.0031.1546	8682241840199	Locking Plate, Hand, Thumb phalanx, 7 holes
152.0031.1547	8682241840205	Locking Plate, Hand, Thumb metacarpal, 6 holes
152.0031.1548	8682241840212	Locking Plate, Hand, grid, 4 holes, Left
152.0031.1549	8682241840229	Locking Plate, Hand, grid, 4 holes, Right
152.0031.1550	8682241840236	Locking Plate, Hand, grid, 6 holes, Left
152.0031.1551	8682241840243	Locking Plate, Hand, grid, 6 holes, Right
152.0031.1552	8682241840250	Locking Plate, Hand, grid, 8 holes, Left
152.0031.1553	8682241840267	Locking Plate, Hand, grid, 8 holes, Right



Mini Screws – Titanium Ti6Al4V

152.2210.2004	8682241829439	Screw, mini, Ti6Al4V, D 2.0 x L 4 mm
152.2210.2005	8682241821488	Screw, mini, Ti6Al4V, D 2.0 x L 5 mm
152.2210.2006	8682241821495	Screw, mini, Ti6Al4V, D 2.0 x L 6 mm
152.2210.2007	8682241821501	Screw, mini, Ti6Al4V, D 2.0 x L 7 mm
152.2210.2008	8682241821518	Screw, mini, Ti6Al4V, D 2.0 x L 8 mm
152.2210.2009	8682241821525	Screw, mini, Ti6Al4V, D 2.0 x L 9 mm
152.2210.2010	8682241821532	Screw, mini, Ti6Al4V, D 2.0 x L 10 mm
152.2210.2012	8682241821549	Screw, mini, Ti6Al4V, D 2.0 x L 12 mm
152.2210.2014	8682241821556	Screw, mini, Ti6Al4V, D 2.0 x L 14 mm
152.2210.2016	8682241821563	Screw, mini, Ti6Al4V, D 2.0 x L 16 mm
152.2210.2018	8682241821570	Screw, mini, Ti6Al4V, D 2.0 x L 18 mm

Mini Locking Screws – Titanium Ti6Al4V

152.2210.2105	8682241821587	Screw, mini Locking, Ti6Al4V, D 2.0 x L 5 mm
152.2210.2106	8682241821594	Screw, mini Locking, Ti6Al4V, D 2.0 x L 6 mm
152.2210.2107	8682241821600	Screw, mini Locking, Ti6Al4V, D 2.0 x L 7 mm
152.2210.2108	8682241821617	Screw, mini Locking, Ti6Al4V, D 2.0 x L 8 mm
152.2210.2109	8682241821624	Screw, mini Locking, Ti6Al4V, D 2.0 x L 9 mm
152.2210.2110	8682241821631	Screw, mini Locking, Ti6Al4V, D 2.0 x L 10 mm

152.2210.2112	8682241821648	Screw, mini Locking, Ti6Al4V, D 2.0 x L 12 mm
152.2210.2114	8682241821655	Screw, mini Locking, Ti6Al4V, D 2.0 x L 14 mm
152.2210.2116	8682241821662	Screw, mini Locking, Ti6Al4V, D 2.0 x L 16 mm
152.2210.2118	8682241821679	Screw, mini Locking, Ti6Al4V, D 2.0 x L 18 mm



	CAT #	Item	Qty
①	IN9817	Long K-Wires Ø1.5 mm x 150 mm	3
②	IN2002	Depth Gauge	1
③	IN2003	Bending Irons	2
④	IN2004	2.4mm Drill Guides (Yellow)	2
⑤	IN2005	3.2mm Drill Guides (Magenta)	2
⑥	IN2006	Soft-Handled Torque Driver	1
⑦	IN2007	Soft Handled Torque Driver for 2.4mm (T10) & 3.2mm (T15) Torx Shafts	1
⑧	IN2008	Torque Adapter	1
⑨	IN2009	Threaded Guide for 2.4mm Drills	2
⑩	IN2010	Threaded Guide for 3.2mm Drills	2
⑪	IN2011	2.4mm Torx Driver Shaft (with screw pick up)	1
⑫	IN2012	3.2mm Torx Driver Shaft (with screw pick up)	1
⑬	IN2013	2.4mm Calibrated Drill with AO Drive	2
⑭	IN2014	3.2mm Calibrated Drill with AO Drive	2
⑮	IN2015	Trial Implant Caddies	3
⑯	IN2016	MIS Drill Jig	1
	Various	Trial Plates – Fore, Mid, Hind	18
	Various	Trial Wedges – Fore	4
	IN2040	Instrument Case with Poly Holders	1

Minimally Invasive System

MI Burrs and Cutters

Burrs and Cutters

CAT #	Item	Size
IN9871	Shannon Isham Burrs	Ø2.0 x 12mm
IN9872	Shannon Standard [STD]	-
IN9873	Shannon Corta	-
IN9874	Conical Burr	Ø3.1 x 13mm
IN9875	Conical Burr	Ø4.1x 13mm
IN9876	Cylindrical Burr	Ø2.0 x 20mm
IN9877	Cylindrical Burr	Ø2.5 x 14mm
IN9878	Cylindrical Burr	Ø3.1 x 22mm
IN9879	Wedge Burr Sterile packed - Single Use	Ø3.1 x 13mm
IN9880	Wedge Burr	Ø4.1 x 13mm
IN9881	Brophy Burr	Ø5.0 x 15mm

Wedge
Burr


Instruments

CAT #	Item	Size
IN9884	Scalpel Handle for Micro Blade	-
IN9885	Double Ended Bone Rasp	-
IN9886	Straight Periosteal Elevator	-
IN9887	Curved Periosteal Elevator	-

Brophy
Burr

Motor Drive Unit

CAT #	Item	Size
IN9891	High Speed Electronic Drive System with Irrigation & Foot Paddle Control	-



Wire-Fix Compression / Distraction Forceps & Pad Distractor Forceps

Wire-Fix Compression / Distraction Forceps



Pad Distractor Forceps



APPENDIX

Ref	Item / Description	Code
A	Osteotomy Instrument Set	IN9854
B	Osteotomy Surgery Pack – TARGET™ Type 1 & Type 2	MD164
C	TARGET™ S3 Drill & K-Wire User Guide	MD163
D	Harpoon™ User Guide	MD169
E	MEGA™ Staple Locking Plate – Custom STG	MD157
F	MEGA™ Compression Staple Locking Plate User Guide	MD162
G	MEGA™ Staple Locking Plate Surgery Pack	MD165
H	Serrated Compression Staple Size Chart	MD161
I	Combined Staple Instrument Set	IN9855
J	Fusion Instrument Set	IN9856
K	1st MTP Fusion D-Tach™ System – Surgery Pack	M167
L	HV Xpanding Spring-Clip™ Instrument Set	IN3026
M	HV Xpanding Spring-Clip™ Surgical Guide	MD170
N	XpandComp™ HF Instrument Set	IN3019
O	XpandComp™ TTC Instrument Set	IN3018
P	Locking Plates and Screws Instrument Set	IN9800
Q	Locking Plates for Fore, Mid, and Hindfoot Indications	MD164
R	Small Extremity Plates and Screw Systems	MD166
S	Wire Fix Distraction and Compression Forceps	IN9882
T	Pad Distraction Forceps	IN9883
U	Optional Sterile Disposable Surgery Sets:	
	• TARGET™ 1 S3	DI1001
	• TARGET™ 2 S3	DI1002
	• MEGA™ Staple	DI1003
	• MTP Plate	DI1004
V	MIS Instrument Set	IN9895
W	MIS Surgical Technique Guide	MD158

MD 160 Rev 5 Dec 2025



71 Jay Avenue, Thornaby-on-Tees
North Riding of Yorkshire TS17 9LZ
United Kingdom.

tel: +44 (0) 1642 769080
email: enquiries@implantsinternational.com
web: www.implantsinternational.com

EC	REP
----	-----

ARAZY GROUP (IRELAND) LTD
19 BAGGOT STREET
LOWER DUBLIN 2 D02
X658 IRELAND PH: +353
21 731 9558

ireland@arazygroup.com www.arazygroup.com